


# RSNA Intracranial Hemorrhage Detection

Keep Digging Gold  
5<sup>th</sup> Place Solution  
November 2019


kaggle




## Team



Tien-Dung Le  
(tarobxl)



Tri Duc Nguyen Tang  
(KeepLearning)



Bac Nguyen



Toan Duc Bui





Anjum Sayed  
(datasaurus)


Preprocessing:  
Creating 3 color  
channels from  
DICOM images

## #1: 3 CT Window Method

Brain window [40, 80]  
Bone window [600, 2800]  
Subdural window [75, 215]



Crop image



Linear window parameters: [width, level]

## #2: Adjacent Slice Method



Reconstruct 3D volume  
from study ID & patient  
position:

*n-1*th image

*n*th image

*n+1*th image

Apply subdural window



Crop image




Image: ID\_0fa314037,  
intraparenchymal

## Overview of models used

Preprocessing Method	Network Architecture	Image Resolution	Score Before Post processing	Score After Post processing
1*	ResNet-18	512x512	0.060	0.054
1*	ResNet-34	512x512	0.058	0.052
1*	ResNet-50	512x512	0.058	0.052
1*	Inception-V3 & Deep Supervision	512x512	0.063	0.053
1	DenseNet-169	512x512	0.055	0.049
1	EfficientNet-B3	300x300	0.055	0.050
2	ResNet-50	512x512	0.054	0.051
2	EfficientNet-B0	224x224	0.054	0.051
2	EfficientNet-B5	456x456	0.048	0.048

\* No cropping on the 3 window preprocessing

**Final private leaderboard score: 0.04561**

## Training parameters

### Cross validation scheme

- 5 fold CV grouped by Patient ID

### Loss

- Weighted binary cross-entropy loss per class (weights = 2 for any and 1 for the other classes)

### Image augmentation

- Horizontal flip, elastic transform, grid distortion, optical distortion, shift-scale-rotate, random resized crop, random rotation


### Optimizers & learning rate:

- AdamW optimizer, initial learning rate: 0.01 or 0.001, weight decay: 0.01
- ReduceLROnPlateau or Cosine annealing, both with early stopping (patience=3)

### Test time augmentation (TTA):

- Identity, horizontal flip, rotate +/- 10 degrees

## Postprocessing




Create a binary classifier for each of the 6 classes using the following features:

- The original prediction  $p_0$
- The predictions of the previous image  $p_{\text{prev}}$  and the next image  $p_{\text{next}}$
- The statistical features for all images before and all the images after the image in the same studyID: number of images, mean/std/skew of the original predictions.

Binary classifiers built using H2O AutoML & OOF predictions

## Model Stacking




### Neural Network:

- Input features:  $54 = 9 \text{ models} \times 6 \text{ classes} \rightarrow \text{reshape } (9, 6, 1)$
- Using CNN:
  - kernel  $(\text{NUM\_MODELS}, 1)$  to learn the correlation between 9 models
  - kernel  $(1, \text{NUM\_CLASSES})$  to learn the correlation between 6 classes
- Mix-up augmentation works pretty well

### Light gradient boosting machine (LightGBM):

- We build 6 separate models for 6 classes

## ROC AUC Performance



Results from validation of  
fold 1 from the  
EfficientNet-B5 model

Classes:

- 0. Any
- 1. Epidural
- 2. Intraparenchymal
- 3. Intraventricular
- 4. Subarachnoid
- 5. Subdural

Similar study achieved average ROC-AUC scores in the region of 0.93. See figure 3d of:  
*Precise diagnosis of intracranial hemorrhage and subtypes using a three-dimensional joint convolutional and recurrent neural network*, Hai et al. April 2019



Important and  
Interesting  
Findings

- Overfitting was a commonly reported issue:
  - Most models only needed approx. 5-7 epochs to reach the best validation score
  - Early stopping was essential
- Both pre-processing methods have strengths & weaknesses:
  - Method #1 benefitted most from the post-processing
  - Method #2 already captured some of the spatial detail but could only utilize one CT window
  - Ideal for model ensembling
- Treating this as a 3D problem gives the best results:
  - Using the metadata effectively was key
  - Worth bearing in mind when labelling future radiological datasets

kaggle™